**B&B DENTAL**  
IMPLANT COMPANY

# PHYSIO BRUSHLESS 3000

## USER MANUAL

**CE 0425**

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ita.

## PHYSIO BRUSHLESS 3000

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## RECOMMENDATIONS

It is very important that this manual is kept with the equipment for future reference. If the equipment is sold or handed over to another user, make sure to hand the manual over to the new owner so that they can learn about the equipment operation and relevant warnings.

These warnings are provided for safety reasons. Read them carefully before installing or operating the equipment.

Warning about environmental protection

Old electrical appliances contain valuable materials, they are not normal household waste! We encourage our customers to help protect the environment and resources and to deliver this equipment to authorised collection centres, if they are available in the area.

**B. & B. Dental S.r.l. - Via San Benedetto, 1837 - 40018 San Pietro in Casale (BO) Italy**

## SAFETY PRECAUTIONS

or

These symbols indicate that the user should refer to the supplied manual.

Type BF (floating applied piece)

For safety reasons, it is essential that the equipment can only be connected to a power outlet with a good earth connection.

Insert contra-angle, drills and any other tool while the motor stopped.

- Use rotating tools (contra-angles, drills or other tools) in perfect working order.
- Use contra-angles suitable for implantology.
- Have the contra-angles checked by competent personnel authorised by the manufacturer on a yearly basis.
- Use the drills specified by your implant supplier.
- Closely observe all instructions given by the manufacturers of the contra-angles and any tools applied to the equipment when selecting maximum speed, maximum torque, clockwise and anticlockwise rotation.
- Connect the equipment to outlets with a solid earth circuit.
- Avoid arranging the connection cables of drill and foot control in tight bends.
- NOTE: To avoid overheating the drill during procedures that require high torque, stop the drill for 10 seconds every 30 seconds of operation.
- Report any serious incidents to the manufacturer and competent authorities.

# STANDARD EQUIPMENT

PHYSIO BRUSHLESS 3000 is made up of the following parts:

- Control unit (Central unit)
- Foot control with cable and connector
- Micromotor with autoclavable cable and connector
- Micromotor protective cap for sterilisation
- Autoclavable micromotor holder
- Hanging rod for saline bag
- Sterile 1-way irrigation tubing
- Power cable


**Accessories not supplied by the manufacturer:**

Contra-angle

Drills

Implants

Consumables: saline solution, disposable irrigation tubings.

Micromotor protective cap for sterilisation

Irrigation tubing

Power cable

Hanging rod for saline bag

## DESCRIPTION

---

Physio Brushless 3000 is a surgical micromotor system for implantology with physio-dispenser featuring electronic speed and torque control, that combines ease of use and high performance.

Both hardware and software were entirely designed by B. & B. Dental S.r.l. with the valuable contribution of RTC Sistemi s.n.c., a well-established name in the industry.

A multi-function foot control allows user to control all functions without touching the front panel.

It uses an autoclavable brushless induction motor with a top speed of 40,000 RPM that accepts handpieces in compliance with ISO 3964.

The peristaltic pump is mounted outside the unit to prevent fluids dripping into the central unit in the event the irrigation tubing is damaged. Its all-metal construction ensures great sturdiness and durability.

The multi-function solid-state foot control is extremely innovative and has no mechanical moving parts, while providing the same functions as the other foot controls. Built on a solid stainless steel base, it sports a stylish oval shape and is fully covered with resilient synthetic rubber that ensures proper sealing against fluid drips. Only 1 centimetre high above the floor, it is equipped with a force sensor that determines motor speed varying it progressively and proportionally to the pressure exerted by user foot.


## INTENDED USE

---

This device is intended for use in dental clinics to perform all necessary operations in implantology procedures. It may be used only by medical personnel specialising in dentistry, qualified and authorised to practice the profession.

The device administers energy to other medical devices (handpiece and surgical drills) and belongs to Class IIa according to Rule 9 of Annex IX to Directive 93/42/EEC. As for electrical safety, the equipment is classified as class I type BF.

# PROGRAMS

PHYSIO BRUSHLESS 3000 offers two program management modes and the choice to use or not the rheostat function of the motor foot control.

## ① PROGRAM

In the 5-program mode, PHYSIO BRUSHLESS 3000 offers the possibility to customise 5 programs by setting the parameters: drill speed and torque, pump flow rate and direction of rotation of the drill.

In the phase mode, PHYSIO BRUSHLESS 3000 offers the possibility to customise 5 programs, with 4 operating phases for each program, by setting the 4 parameters - speed, drill torque, pump flow rate and direction of rotation of the drill - in each operating phase. This way, users can set 5 different presets for implants with different parameters.

For example: use the PROGRAM button to select the number for the implant to be set. PHASE1 will be selected automatically.

PHASE1 - select the 4 operating parameters to mark the implant site.

Press the button no. 2 of the foot control (see foot control description) to move to PHASE2 that lets you select the 4 operating parameters for pilot milling.

Press the button no. 2 of the foot control (see foot control description) to move to PHASE3 that lets you select the 4 operating parameters for milling and tapping.  
Press the button no. 2 of the foot control (see foot control description) to move to PHASE4 that lets you select the 4 operating parameters for screwing.

### ② **RPM (drill revolutions per minute)**

Motor speed can be adjusted using the appropriate buttons, in 22 steps in a range from 20g/min to 2,000g/min, with a tolerance of  $\pm 10\%$ .

### ③ **TORQUE**

Controlled torque can be selected. There are 24 settings available and the value is shown on the display in Ncm, with a tolerance of  $\pm 10\%$ . For example, at the output of the contra-angle 20:1, you can select torque values of 3Ncm, 6.9,.....69 Ncm in 3Ncm steps, using the appropriate buttons. In the MAX position, there is not torque limitation.

### ④ **CALIBRATION**

With the contra-angle inserted in the motor and placed on its holder, press this button to calibrate the contra-angle in all available speeds.

### ⑤ **PUMP**

The PUMP button lets you choose the flow rate at which saline solution is delivered among 3 values: min, mid, max, or OFF to stop delivery.

### ⑥ **DRILL**

Starts and stops the drill.

### ⑦ **REVERSE**

Press to reverse the rotation of the drill. An acoustic signal occurs when REV mode is on.

### ⑧ **REDUCT**

Available reduction ratios for the contra-angles are as follows: 20:1, 20:1Wh.

The contra-angle 20:1Wh is the W&H precision contra-angle mod.Wi-75 R/KM or equivalent.

Torque calibration is carried out with the above-mentioned W&H contra-angle that provides the best torque control accuracy. When you use lower quality contra-angles with a 20:1 ratio, select 20:1 as reducer contra-angle.

#### **Torque and RPM**

In order to achieve the utmost accuracy in limiting the torque delivered by the motor in the lower torque range, motor RPM will be limited as outlined below:

> with Reduct 20:1

> with torques from 3Ncm up to 12Ncm maximum motor speed is 1250 RPM

> with torques from 15Ncm up to 36Ncm maximum motor speed is 1750 RPM

> with torques higher than 39Ncm maximum motor speed is 2000 RPM.

## ⑨ CONFIGURATION

Pressing the Config button on the display shows the configuration menu. Pressing the Program button in this menu lets you choose among 4 combinations of phase and rheostat modes.

Once you have obtained the desired combination, as shown on the display, press the Config key to store your selection and exit the configuration menu.

## FOOT CONTROL

1. The left button selects pump flow rate.

2. The rear central button (under the handle), in the 5-program mode, selects the program to run.

In the phase mode, it selects the PHASE to run, or set, (PHASE1,2,3,4) within the program you will have selected using the PROGRAM button on the central unit panel.

3. The right button reverses drill rotation.

4. The central front and oval button operates the drill and, if selected, the pump. If the RHEOSTAT function is in the OFF state, it allows the simple start or stop of the drill at the speed indicated on the display; instead, if the RHEOSTAT function is activated, therefore in the ON state, the force sensor allows you to operate the drill with variable speed by exerting more or less pressure on the corresponding button, where the maximum achievable speed remains the one set and displayed.

## INSTALLATION AND SET UP

- Remove all equipment components from the packaging without damaging it and keep the package in case you need to transport the equipment in the future.
- Check the equipment for transport damage. Report any damage to your supplier within 24 hours of receipt.
- Check that the mains voltage selector switch is set to the voltage that the equipment will be running on.
- Connect foot control and drill to the central unit.
- Insert the contra-angle to be used on the motor and place it on the supplied holder.
  - The motor accepts any contra-angle compliant with ISO 3964 standard
  - Recommended reduction ratio is 20:1.
- Connect the power cable and switch on the equipment.

• Press the CALIBRATION button and wait until calibration is completed, as indicated by the display.

**Example of "5-program mode"**

- From PROGRAM 1, select the REDUCT setting corresponding to the contra-angle inserted; set RPM, TORQUE and PUMP to the values for the first of your work phases.
- Repeat the above steps in the following PROGRAMS for the next work steps as needed.
- When finished with the above procedures, the parameters chosen for each program will have been stored automatically. Press the PROGRAM button to review and edit them as required.
- If the equipment is switched off and back on again, program parameters will be those set upon last access.

## REAR PANEL AND PERISTALTIC PUMP

The following figures show the rear panel with the clip closed (Fig.1) and with the clip open (Fig.2)

Fig.1


Fig. 2

**Please note: below is the key to some label information:**

- Fuse: 1.6 A fuses for 115 V; 1 A fuses for 230 V
- Intermittent operation Ton/Toff: 30 sec ON and 10 sec OFF
- Manufactured by: Manufacturer's details.

## INSERTING THE IRRIGATION TUBING

Starting condition - clip released - is shown in fig.3. Route the irrigation tubing as shown in fig.4, and then fig.5, with fig.6 showing the final condition with the clip closed.

Fig.3


Fig. 4

Fig.5


Fig. 6

We recommend using the "INFUSET" irrigation tubing supplied by B. & B. Dental. However, the device will accept any silicone tubing with an internal diameter of 3.5 mm, an external diameter of 5.0 mm and a length between 10.5 and 12.0 cm (see fig. 7)


Fig.7

# MAINTENANCE

## CLEANING AND STERILISATION

- All Physio Brushless 3000 components can be cleaned with a cloth moderately wetted with alcohol or moistened with water. The screen-printed graphics of the central unit panel should only be cleaned with a cloth dampened with water only.

- The drill, along with cable and connector, and the contra-angle holder are autoclavable at a maximum temperature of 135°C.

The drill can undergo a minimum of 1000 sterilisation cycles in an autoclave.

Autoclaving temperature: 134°C for at least 3-4 minutes; Pressure: 2.3 bar; Relative humidity: 100%; Total cycle time: 20 min.

## IMPORTANT INSTRUCTIONS FOR MICROMOTOR STERILISATION

The life of the micromotor is affected by the number of sterilisation cycles. As a result, the following precautions must be observed or the warranty will become null and void.

**Always insert the supplied cap on the micromotor, wrap the cable without creating bends and insert it in the appropriate bag before placing it in the autoclave. Use only distilled water in the autoclave, do not use demineralised water.**

## ROUTINE MAINTENANCE

For the maintenance of the contra-angle, observe the manufacturer's instructions.

- At least every 6 months, check power cable and drill cable for abrasion or breakage. If you find any such damage, have the component replaced at an approved centre.

- In the event of faults or malfunctions during operation, switch off the equipment immediately and contact your dealer. Report the type of defect and locate the authorised service centre in charge of your area. If there is none, contact B. & B. Dental directly.

- Under the regulations in force, the equipment must undergo a periodic safety inspection every 2 years.

# TECHNICAL SPECIFICATIONS

Manufacturer	B. & B. Dental S.r.l. Via San Benedetto, 1837 - 40018 San Pietro in Casale (BO) - Italy Tel. +39 (0) 51.81.13.75 - Fax +39 (0) 51.666.94.00 www.bebdental.it
Model	PHYSIO BRUSHLESS 3000
Dimensions of central unit	255 x 255 x 110 mm
Central unit	ABS and aluminium
Drill	Aluminium
Foot control	Stainless steel and rubber
Weight	6 kg (2.6 Kg central unit)
Micromotor	91,2 W - Nominal Torque 18,7 mNm - 44,300 RPM - Sterilisable
Noise	<65 dBA
Tolerance in indicated speeds	±10%
Tolerance in indicated torques	± 10%
Supply voltage	AC 115 / 230 V ± 10%      50/60 Hz
Classification according EN 60601-1	Class I type BF
Classification according Directive 93/42/ECC	Medical device: IIa , Rule 9
Use	Intermittent: 30 sec ON and 10 sec OFF
Environmental operating conditions	Temperature: +10°C to +40°C NOT SUITABLE for use in the presence of flammable anaesthetics and oxygen
Storage and transport	From 5°C to 70°C
Rated power	100 VA
Fuses	5 x 20 mm
Power supply 115V	2 x 1.6 HT
Power supply 230V	2 x 1 HT

## TROUBLESHOOTING

THE CENTRAL UNIT WILL NOT TURN ON	<ul style="list-style-type: none"> <li>• Check the connection between unit and mains.</li> <li>• Check the condition of the fuses.</li> </ul>
THE FOOT CONTROL WILL NOT WORK	Check the connection between foot control and unit, if the connection is correct have the foot control serviced.
THE DISPLAY SHOWS UNUSUAL CHARACTERS	If the display shows unusual fonts, the user can proceed with the operations ongoing, since the device continues to work using the settings previously displayed. In alternative, the user can turn off and on the device, so that the display still shows the previous settings, or the user can perform the system reset by switching the device off and on, while pressing the PROGRAM button. In this last case, for a few seconds the display will show the Firmware/ software version and it will return to the factory setting.
THE DRILL WILL NOT RUN	<ul style="list-style-type: none"> <li>• Remove the contra-angle and operate the drill again.</li> <li>• If the drill works properly, have the contra-angle serviced.</li> <li>• If it still will not run, have the entire system serviced.</li> </ul>
THE PERISTALTIC PUMP WILL NOT RUN	<ul style="list-style-type: none"> <li>• Ensure that the pump rotor is running and the irrigation tubing is correctly in place.</li> </ul>

## WARRANTY TERMS AND CONDITIONS

- The product is guaranteed for one year from date of invoice.
- The drill motor is guaranteed for six months.
- The warranty does not cover the costs of transport and / or service calls at the customer's premises.
- In the event the equipment is serviced or repaired by unauthorised personnel, the warranty will become null and void.
- The purchaser has no right to compensation for equipment downtime.
- The warranty does not cover any faults that:
  - are caused by wear and tear due to operation, as may be the case with motor brushes and bearings, foot control and drill cables, power cable;
  - are caused by fire or weather events such as lightning, etc.

**To obtain service, please contact us at:**  
**B. & B. Dental S.r.l. - Via San Benedetto, 1837**  
**40018 San Pietro in Casale (BO) - Italy**  
**Tel. +39 (0) 51.81.13.75 - Fax +39 (0) 51.666.94.00**  
**www.bebdental.it**

# DECLARATION OF CONFORMITY

MANUFACTURER:	<b>B. &amp; B. Dental S.r.l.</b> Via San Benedetto, 1837 40018 S. Pietro in Casale (BO) - Italy
MEDICAL DEVICE:	<b>PHYSIO BRUSHLESS 3000</b>
PRODUCT CODE (REF):	<b>PHYSIO BRUSHLESS 3000</b>
CLASSIFICATION ANNEX IX MDD:	<b>CLASS IIa RULE 9</b>
CERTIFICATION ROUTE:	<b>ANNEX V</b>

B. & B. Dental S.r.l. herewith certifies that the surgical micromotor for implantology, marketed with:

Serial number:

It complies with the essential requirements set out in Annex I of Directive 93/42/EEC and subsequent amendments and additions (Directive 2007/47/EC) as specified in Annex II of Directive 93/42/EEC and subsequent amendments and additions (Directive 2007/47/EC) and transitional provisions according to Art. 120 (3) of Regulation (EU) 2017/745.

It also declares that the device complies with the following harmonised standards:

EN 60601-1: Medical electrical equipment. Part 1: General requirements for safety"  
EN 60601-1-2: "Medical electrical equipment. Part 1-2: General requirements for safety -  
Collateral standard: Electromagnetic compatibility. Requirements and tests".

## Harmonized European applicable rules

The list of applicable standards is reported on the Technical File for the product in issue.

Notified Body: <b>ICIM</b> Piazza Don Enrico Mapelli 75 20099 Sesto San Giovanni (MI)	Notified Body N°: <b>0425</b>
Certificate N° : 0425-MED-004030-00	Date of Certificate: <b>17-07-2020</b>
Expiry Date: <b>26-05-2024</b>	Declaration Expiry: <b>26-05-2024</b>

San Pietro in Casale, \_\_/\_\_/\_\_\_\_

Legal Representative  
BIANCOSPINO ANNA

# ACCOMPANYING DOCUMENTS

Below are the EMC emission and immunity tables provided as guidance for the manufacturer, compiled in accordance with article 6 of EN 60601-1-2 based on the tests and results obtained. The manufacturer must report this information in the documentation that comes with the product.

Emission aspects		
Emission tests	Conformity	Electromagnetic environment - Guidelines
RF emissions Cispr 11	Group 1	The product Physio Brushless 3000 uses RF energy only for its internal operation. As a result, its RF emissions are very low and are unlikely to cause interference in nearby electronic devices.
RF emissions Cispr 11	Class A	The product Physio Brushless 3000 is suitable for use in all buildings including household buildings and those directly connected to a low voltage power supply network supplying household buildings. <b>Warning:</b> the equipment is intended for use by qualified personnel only. In addition, this equipment may cause radio interference or cause nearby equipment to malfunction. It may be necessary to implement appropriate measures to reduce such interference, such as changing the orientation or location of the equipment and its cables, or even shielding the site.
Harmonic emissions IEC 61000-3-2	Class A	The equipment can be used in all buildings including household buildings and those directly connected to the public low voltage power supply network supplying household buildings.
Voltage fluctuation/ flicker emissions IEC 61000-3-3	Compliant	

<b>Immunity aspects</b>			
Physio Brushless 3000 has been designed to operate in the electromagnetic environment specified below. The customer or user should ensure that it is used in such an environment.			
<b>Immunity testing</b>	<b>EN 60601-1-2 testing level</b>	<b>Compliance level</b>	<b>Electromagnetic environment - Guidelines</b>
Electrostatic discharge (ESD) EN 61000-4-2	± 8kV when in contact ± 15kV in the air	± 8kV when in contact ± 15kV in the air	Floors must be made of wood, concrete or ceramic. If floors are covered with synthetic material, relative humidity should be at least 30 %.
Electrical fast transients/bursts EN 61000-4-4	±2kV power supply lines	±2kV power supply lines	The quality of the mains voltage should be that of a typical commercial or hospital environment.
Impulses EN 61000-4-5	± 1kV differential mode ± 2kV common mode	± 1kV differential mode ± 2kV common mode	The quality of the mains voltage should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions, and voltage variations on input lines EN 61000-4-11	< 5% $U_T$ (>95% dip in $U_T$ ) for 0.5 cycles 70% $U_T$ (30% dip in $U_T$ ) for 25 cycles < 5% $U_T$ (>95% dip in $U_T$ ) for 5 seconds	< 5% $U_T$ (>95% dip in $U_T$ ) for 0.5 cycles 70% $U_T$ (30% dip in $U_T$ ) for 25 cycles < 5% $U_T$ (>95% dip in $U_T$ ) for 5 seconds	The quality of the mains voltage should be that of a typical commercial or hospital environment. If the user requires continuous operation even during a mains outage, we recommend using an uninterruptible power supply (UPS) unit or batteries to supply the equipment.
Power frequency magnetic field EN 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a commercial or hospital environment.

**RF immunity aspects**

Physio Brushless 3000 has been designed to operate in the electromagnetic environment specified below. The customer or user should ensure that it is used in such an environment

Immunity testing	EN 60601-1-2 testing level	Compliance level	Electromagnetic environment - Guidelines
Conducted RF			
EN 61000-4-6	3 Vrms 150kHz to 80MHz	3 Vrms 150kHz to 80MHz	Portable and mobile RF communication equipment should not be used closer to any part of the equipment, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distances: $d = 1.2 \cdot \sqrt{P}$ 150kHz to 80MHz $d = 1.2 \cdot \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \cdot \sqrt{P}$ 800 MHz to 2.5 GHz $d = 2.3 \cdot \sqrt{P}$ 800 MHz to 2.5 GHz Where P is the maximum output rated power of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).
Radiated RF			
EN 61000-4-3	3 V/m 80MHz to 2.7GHz	3 V/m 80MHz to 2.7GHz	

The field strength from fixed RF transmitters, as determined in an electromagnetic site survey, may be less than the compliance level in each frequency range.

Interference may occur in the vicinity of equipment marked with the following symbol:

Recommended separation distance between portable and mobile radio communications equipment and the equipment.

Physio Brushless 3000 has been designed to operate in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the operator of the equipment can help to prevent electromagnetic interference by ensuring a minimum distance between mobile and portable RF communication equipment (transmitters) and the equipment, as recommended below, in relation to the maximum output power of radio communication equipment.

Rated maximum output power of transmitter (W)	Separation distance at transmitter frequency (m)		
	150kHz to 80MHz $d = 1.2 \cdot \sqrt{P}$	80MHz to 800MHz $d = 1.2 \cdot \sqrt{P}$	800MHz to 2.5GHz $d = 2.3 \cdot \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters with maximum output rated power not listed above, the recommended separation distance  $d$  in metres (m) can be calculated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output rated power of the transmitter in watts (W) according to the transmitter manufacturer.

**Note:**

- (1) At 80 MHz and 800 MHz, the higher frequency range applies.
- (2) These guidelines may not apply in all situations. Electromagnetic propagation is influenced by the absorption and reflection from structures, objects and people.

---

# CONTACTS

## REGISTERED OFFICE & OPERATIONAL 1:

Via S. Benedetto, 1837

40018 S. Pietro in Casale (BO) Italy

## OPERATIONAL HEADQUARTERS 2:

Via Due Ponti, 19

40050 Argelato (BO) Italy

Tel. +39 (0) 51.81.13.75

Fax +39 (0) 51.666.94.00

info@bebdental.it

www.bebdental.it